

## Special Investment for Lithium Disilicate Ceramics.

This investment which was especially developed for processing of Lithium Disilicate Glass Ceramics guarantees...

- perfect surfaces
- a minimal reaction layer after pressing
- high detail reproduction of the pressed objects
- consistent expansion values due to strict quality controls



### CONCEPT VEST™ Powder

Special Investment  
for Lithium Disilicate Ceramics.



25 x 100 g

Item-No. 898/2066

50 x 100 g

Item-No. 898/2065

### Mixing Liquid



1 Liter

Item-No. 898/204

### Zubler® Disposable Plungers

For Pressing of Dental Ceramics  
and **Lithium Disilicate**.

Their excellent quality and exact dimensions guarantee  
the repeatable success of pressings in rings of all sizes.



ø 12 mm **Item-No. 20/0235**

ø 13 mm **Item-No. 20/0238**

ø 16 mm **Item-No. 20/0237**

### Zubler® Flexring System for Lithium Disilicate Ingots

The ribbed surface of the rings ensures faster heat absorption and re-  
lease. Sector lines make it easier to find the pressed object, and damage  
during the divesting process is virtually impossible. The flexring is easy  
to clean and has a long life span.



Flexring set 100 g **Item-No. 20/0228**

Flexring set 200 g **Item-No. 20/0229**

## VARIO BALANCE

The dosage system for the demanding dental lab  
- 10 x higher accuracy compared to dosing  
with a graduated cylinder.

- 200 customizable programs
- fast and easy to the highest precision
- Extremely consistent results
- Highly accurate, maximum error rate (< 1%)



Item-No. 898/420